

On Customer Site

## Visit Report

**Reference :**

|                         |                  |                    |                      |
|-------------------------|------------------|--------------------|----------------------|
| Commande BSTe France :  |                  | Call : E1275       |                      |
| Customer :              | Valetiq          | System :           | TubeScan DS 430 4000 |
| Customer Contact :      | M. Arthur Caniau | SN :               | 42.2234 / ASO087513  |
| Customer Order number : | 222143           | Machine Nickname : | Omega 3              |

**Visit Purpose :**

|                        |                                         |                |                                     |                 |                                     |
|------------------------|-----------------------------------------|----------------|-------------------------------------|-----------------|-------------------------------------|
| Mechanical Preparation | <input type="checkbox"/>                | Installation   | <input checked="" type="checkbox"/> | Start-up        | <input checked="" type="checkbox"/> |
| Maintenance            | <input type="checkbox"/>                | Under Contract | <input type="checkbox"/>            | Out of Contract | <input type="checkbox"/>            |
| Troubleshoot           | <input type="checkbox"/>                | Under Warranty | <input type="checkbox"/>            | Out of Warranty | <input type="checkbox"/>            |
| Instruction            | <input checked="" type="checkbox"/> Yes |                |                                     |                 |                                     |
|                        | <input type="checkbox"/> No             |                |                                     |                 |                                     |

|             |            |  |  |  |  |
|-------------|------------|--|--|--|--|
|             | 02/08/2022 |  |  |  |  |
| Travel To   | 3.5        |  |  |  |  |
| Start Time  | 15h00      |  |  |  |  |
| Lunch Time  | 0.5        |  |  |  |  |
| End Time    | 18h00      |  |  |  |  |
| Travel From | 0.5        |  |  |  |  |

|         |                                             |
|---------|---------------------------------------------|
| Details | Installation et mise en service du TubeScan |
|---------|---------------------------------------------|

Purpose of the visit completed :      Yes ☐      No ☒

Need a new visit :      Yes ☒      No ☐

*On Customer Site*

---

## Customer Details for NYQUIST Systems

---

Customer name and place : [Valetiq, rue Berthelot, 59860 Bruay-sur-l'Escaut – France](#)

Contact person at customer : [M. Arthur Caniau \(Responsable Technique\)](#)

Customer Mail & Mobile/Landline : [Arthur.caniau@valetiq.com](mailto:Arthur.caniau@valetiq.com) / +33 3 27 30 97 92 / +33 6 88 82 60 86

System Type : [TubeScan DS 430 4000/X104iC-0.1111](#)

SN : [42.2234](#)


Machine Manufacturer/Model/Nickname : [Omega 3](#)

Date : [03/08/2022](#)

Name : [Richard Relin](#)

**Signatures :**

BST eltromat FSE



Richard Relin 2022-08-04T13:01:49+02:00

Customer



Arthur Caniau 2022-08-04T12:45:29+02:00

On Customer Site

---

## Technical datas

---

Computer name : [Tubescan422234](#)

Total Physical Memory : [16 GB](#)

Free Space in C Drive (System) : [144 GB](#)

OS : [Win 10 Pro](#)

Screen Resolution Touch : [1024x768](#)

Screen Resolution Full View : [2560x1440](#)

Computer Model : [Optiplex 7090](#)

Date : [2022](#)

Version : [7.4.16.36](#)

### **Remote Maintenance (Teamviewer)**

IP address :

Subnet mask :

Default gateway :

Preferred DNS server :

Alternate DNS server :

**Teamviewer ID / Version :** [561 004 738/ v12](#)

*On Customer Site*

---

Notes

---